Instructions for Use

Hrvatski/Croatian **p2**

Cesky/Czech **p3**

Dansk/Danish **p4**

Nederlands/Dutch **p5**

English **p6**

Eesti/Estonian **p7**

Suomi/Finnish **p8**

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Deutsch/German **p10**

Elinika/Greek **p11**

Magyar/Hungarian **p12**

Íslenska/Icelandic **p13**

Italiano/Italian **p14**

한국인 /Korean **p15**

Lietuviška/Lithuanian **p16**

Norsk/Norwegian **p17**

Polski/Polish **p18**

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Română/Romanian **p20**

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INSTRUCTIONS FOR USE:**BUILD-IT® LIGHT CURE CORE MATERIAL INSTRUCTIONS FOR USE****Indications:**

- Core build material on vital teeth
- Core build up material on non-vital teeth
- Restorative material

Contraindications:

- Patients that may be sensitive to any of the ingredients

INSTRUCTIONS FOR USE:

1. Isolate the tooth with a rubber dam or equivalent isolation.
2. Remove all caries according to standard techniques.
3. For deep cavities on vital teeth cover the floor adjacent to the pulp with calcium hydroxide or glass ionomer.
4. Apply a bonding agent to the prepared dentin and enamel according to manufacturer's instructions. Refer to the instructions of your preferred bonding agent. Bond1® SF Solvent-Free SE Adhesive, NanoBond® Adhesive System, Bond-1® Primer Adhesive, and Bond-It® Enamel/Dentin Bonding System are recommended.
5. Apply and adapt *Build-It* Light Cure Core Material to the preparation and post, if present, and form the appropriate anatomy with an instrument. It is recommended that the material be tamped into the preparation with an instrument to eliminate any potential voids prior to curing.
6. Light cure *Build-It* Light Cure Core Material for 20 seconds per surface with a curing light with an output in excess of 600mW/cm². *Build-It* Light Cure Core Material has a direct depth of cure up to 10mm. Ensure that the light tip is in close proximity to the material to achieve maximum curing. After the first 5 seconds of curing of the occlusal surface, place the light tip in direct contact with the core material to maximize depth of cure.
7. Prepare and finish the core using diamond and carbide burs. Ensure that a ferrule of 2mm is remaining on natural tooth structure.

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