<u>Hrvatski/Croatian</u>	p2
<u>Česky/Czech</u>	p3
<u>Dansk/Danish</u>	p4
<u>Nederlands/Dutch</u>	p5
<u>English</u>	p6
<u>Eesti/Estonian</u>	p7
<u>Suomi/Finnish</u>	p8
<u>Français/French</u>	p9
<u>Deutsch/German</u>	p10
<u>Elinika/Greek</u>	p11
<u>Magyar/Hungarian</u>	p12
<u>Íslenska/Icelandic</u>	p13
<u>Italiano/Italian</u>	p14
<u>Lietuviška/Lithuanian</u>	p15
<u>Norsk/Norwegian</u>	p16
<u>Polski/Polish</u>	p17
<u>Português/Portuguese</u>	p18
<u>Română/Romanian</u>	p19
<u>Srpski/Serbian</u>	p20
<u>Slovenčina/Slovakian</u>	p21
<u>Slovenščina/Slovenian</u>	p22
<u>Español/Spanish</u>	p23
<u>Swenska/Swedish</u>	p24
<u>Türkçe/Turkish</u>	p25

UPUTE ZA UPORABU

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INSTRUCTIONS FOR USE:

In order to assist you when using the FibreKor™ Post System, we have included the following **Helpful Hints:**

- **To enhance the resistance fracture of the core and restoration, 2mm of tooth structure should be left above the gumline prior to and after preparation.**
- *FibreKor* Posts may be sterilized. Use cold sterilants such as glutaraldehyde.
- *FibreKor* Drills must be sterilized.
- *FibreKor* Posts are easily cut with a diamond separating disk.
- Cement may be carried to the canal via a Lentulo Spiral instrument to ensure good coverage of cement in canal.
- Finish line or margin of preparation should remain on solid tooth structure, not on core build-up.
- Resin cements and core materials provide the best result when used with the *FibreKor* Post System.
- If retreatment of the apex is required, the *FibreKor* Post can be removed by drilling into it with the corresponding size of *FibreKor* Drill.

FibreKor Post Drill	ISO File	Gates-Glidden	Peeso
1.00mm	#70 Red (1.02mm)	#4 (1.135mm)	#3 (1.02mm)
1.125mm	#90 White (1.20mm)	#4 (1.135mm)	#4 (1.23mm)
1.25mm	#100 Yellow (1.32mm)	#5 (1.226mm)	#4 (1.23mm)
1.375mm	#100 Yellow (1.32mm)	#6 (1.446mm)	#5 (1.50mm)
1.50mm	#120 Blue (1.52mm)	#6 (1.446mm)	#6 (1.58mm)
1.75mm	#140 Black (1.70mm)	#6 (1.446mm)	#6 (1.70mm)

Instructions for Placing The *FibreKor* Post

1. Place an appropriately sized *FibreKor* Drill against a periapical radiograph to determine the required width of the *FibreKor* Post.
2. Using this same drill, determine the depth of the post space preparation, leaving at least 4mm of gutta percha at the apical end. Place rubber endodontic stops on the drills corresponding to the desired depth to avoid over-preparing the post space.
3. Using either Gates-Glidden drills or Peeso Reamers, remove the gutta percha from the canal. Leave at least 4mm of gutta percha at the apical end.
4. Start the post space preparation (the canal) by first using the smallest drill and then successively larger drills until the desired width is achieved. In the case of the 1mm *FibreKor* Post, a #2 Peeso Reamer or a #4 Gates-Glidden Drill may be used prior to the corresponding *FibreKor* Drill. (see chart)
5. Rinse out the canal to remove the debris and dry it with your air syringe and paper points.
6. Trial seat a *FibreKor* Post in the canal that corresponds to the size of the last drill used. A radiograph may be taken at this point to verify that the post is fully seated.
7. Remove the post and, using a water-cooled high-speed diamond instrument, trim the post to the desired occlusal height.
8. Clean the post with an alcohol wipe and dry it with a blast of oil-free air.
9. Acid-etch the canal with 37% phosphoric acid for 20 seconds.
10. Rinse the canal thoroughly and dry using paper points. Leave the dentin slightly moist.
11. Introduce Bond-1™ Primer/Adhesive (or other suitable bonding agent) into the canal using a syringe or a microbrush. If a microbrush is used, two or more applications may be required. Gently dry the canal with a stream of oil-free air and use a paper point to remove excess resin from the canal. Do not light cure at this point.
12. Mix equal amounts of Cement-It™ Universal C&B™ Cement, Lute-It™ Dual Cure Luting Cement base and catalyst, or Build-It™ FR™ Core Build-up Material (or other suitable dual cure resin cement). Coat the walls of the canal with the cement using a Lentulo Spiral.
13. Coat the Post with *Cement-It Universal C&B* Cement, *Lute-It* Cement or *Build-It FR* Core Build-up Material and insert it into the canal with gentle finger pressure, allowing the excess cement to vent.
14. Remove the excess cement with a plastic instrument.
15. Initiate the set of the dual cure cement by light curing the top of the canal for 10 seconds. *Cement-It Universal C&B* Cement, *Lute-It* Cement, or *Build-It FR* Core Build-up Material will fully set in about 4 minutes.
16. Using a matrix band or core form, a core can be fabricated by bonding a composite core material (*Build-It FR* Core Build-up Material) to the post and remaining dentin.
17. Shape the cured core material to receive the final restoration using diamond and carbide burs.

Bond-1, Cement-It, Lute-It, Build-It and Pentron are registered trademarks of Pentron™ Clinical Technologies, LLC. Universal C&B and FR. trademarks of Pentron Clinical Technologies, LLC.

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