**Instructions for Use**

<b><u>Hrvatski/Croatian</u></b> . . . . .	<b>p2</b>
<b><u>Cesky/Czech</u></b> . . . . .	<b>p3</b>
<b><u>Dansk/Danish</u></b> . . . . .	<b>p4</b>
<b><u>Nederlands/Dutch</u></b> . . . . .	<b>p5</b>
<b><u>English</u></b> . . . . .	<b>p6</b>
<b><u>Eesti/Estonian</u></b> . . . . .	<b>p7</b>
<b><u>Suomi/Finnish</u></b> . . . . .	<b>p8</b>
<b><u>Français/French</u></b> . . . . .	<b>p9</b>
<b><u>Deutsch/German</u></b> . . . . .	<b>p10</b>
<b><u>Elinika/Greek</u></b> . . . . .	<b>p11</b>
<b><u>Magyar/Hungarian</u></b> . . . . .	<b>p12</b>
<b><u>Íslenska/Icelandic</u></b> . . . . .	<b>p13</b>
<b><u>Italiano/Italian</u></b> . . . . .	<b>p14</b>
<b><u>한국인 /Korean</u></b> . . . . .	<b>p15</b>
<b><u>Lietuviška/Lithuanian</u></b> . . . . .	<b>p16</b>
<b><u>Norsk/Norwegian</u></b> . . . . .	<b>p17</b>
<b><u>Polski/Polish</u></b> . . . . .	<b>p18</b>
<b><u>Português/Portuguese</u></b> . . . . .	<b>p19</b>
<b><u>Română/Romanian</u></b> . . . . .	<b>p20</b>
<b><u>Русский язык/Russian</u></b> . . . . .	<b>p21</b>
<b><u>Srpski/Serbian</u></b> . . . . .	<b>p22</b>
<b><u>Slovenčina/Slovakian</u></b> . . . . .	<b>p23</b>
<b><u>Slovenščina/Slovenian</u></b> . . . . .	<b>p24</b>
<b><u>Español/Spanish</u></b> . . . . .	<b>p25</b>
<b><u>Swenska/Swedish</u></b> . . . . .	<b>p26</b>
<b><u>Türkçe/Turkish</u></b> . . . . .	<b>p27</b>

## UPUTE ZA UPORABU

### TempSpan® CMT privremeni cement TempSpan® prozirni privremeni cement

TempSpan privremeni cement na bazi smole, cement bez eugenola namijenjen za održavanje privremenih restoracija. Oba cementa karakterizira naša prikladna dvocijevna šprica sa samomiješajućim vršcima koji vam omogućuju da postavite cement direktno na privremenu restoraciju.

#### Upute za uporabu

##### Indikacije

- Cementiranje privremenih krunica, mostova i nadomjestaka

##### Kontraindikacije

- Permanentno cementiranje konačnih restoracija
- Privremeno cementiranje privremenih navlaka

Na idealnim uvjetima miješanja od 70-77 °F (21-25 °C) i 40-60% relativne vlažnosti, dobiveno je vrijeme rada i vrijeme podešenja s početka miješanja.

##### Privremeni cement

TempSpan® CMT privremeni cement  
TempSpan® prozirni privremeni cement

##### Vrijeme rada

1 - 1½ minute  
≥1½ minute

##### Vrijeme podešenja

2½ - 3 minute  
≤5½ minute

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**TempSpan® CMT Temporary Cement**  
**TempSpan® Transparent Temporary Cement**

TempSpan Temporary Cements are resin based, eugenol-free cements designed for retaining provisional restorations. Both cements feature our convenient double barrel syringe with auto-mix tips that enable you to place cement directly into the provisional restoration.

**Instructions for Use**

**Indications**

- Cementation of temporary crowns, bridges, inlays, and onlays

**Contraindications**

- Permanent cementation of definitive restorations
- Temporary cementation of provisional veneers

At the ideal mixing conditions of 70-77 °F (21-25 °C) and 40-60% relative humidity, the following work time and setting time from start of mixing are obtained.

<b>Temporary Cement</b>	<b>Work Time</b>	<b>Set Time</b>
TempSpan® CMT Temporary Cement	1 - 1½ minutes	2½ - 3 minutes
TempSpan® Transparent Temporary Cement	≥1½ minutes	≤5½ minutes

1. Prior to placing the mixing tip onto the syringe express cement until both sides flow evenly.
2. Dry the prepared teeth and the internal surface of the restoration.
3. Using slight pressure, dispense a thin layer of TempSpan® Temporary Cement directly into the provisional restoration.
4. Within the specified work time, carefully, but firmly, seat the provisional restoration on the preparation site.
5. Instruct the patient to gently occlude into the provisional until the restoration is fully seated.
6. Clean any excess TempSpan® Temporary Cement with a cotton swab prior to setting.
7. TempSpan® CMT Temporary Cement will set in approximately 2½ - 3 minutes from the start of mix. TempSpan® Transparent Temporary Cement will set in approximately 5½ minutes from the start of mix. TempSpan® Transparent Temporary Cement can also be light cured 20 seconds per surface to accelerate the set of the cement.
8. Remove set cement with a scaler or other suitable instrument.

**CAUTION**

TempSpan® Temporary Cements contains (meth) acrylic resins. Uncured (meth) acrylic resins may cause contact dermatitis and damage the pulp. Avoid prolonged contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.

**STORAGE**

Store at ambient temperature.

**SHELF LIFE**

Refer to outer package for expiration date. Do not use after expiration date.

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