

# Correct VPS<sup>®</sup> and Correct Quick<sup>®</sup>

## Vinyl Polysiloxane Impression Materials



The *Correct* System is a full line of vinyl polysiloxane impression materials available in several viscosities, two delivery systems (50ml cartridges and hand mix tubes<sup>†</sup>), and two different set times to satisfy every technique. *Correct VPS* Material offers a 1½–2 minute working time with a 3–4 minute set time, while *Correct Quick* offers a 1 minute working time with a 2 minute set time.

The *Correct* System has a proven track record of providing consistent, reliable results for over 10 years, saving your practice both time and money!

## Accurate, highly detailed impressions!

**Correct VPS/Quick Light Extra Body (~40)\*** – a lower viscosity than traditional light body materials for the ultimate detail.

**Correct VPS/Quick Light Body (~40)\*** – a traditional light viscosity wash material characterized by superior flowability.

**Correct VPS/Quick Medium Body (~45)\*** – a medium viscosity material that produces highly detailed, extremely accurate monophasic impressions; may be used as either a tray material or a light body. Also an ideal, low cost substitute for alginates!

**Correct VPS/Quick Heavy Body (~55)\*** – a high viscosity material that offers excellent hydraulics to drive the wash material into the sulcus.

**Correct VPS Bite™/Quick Bite Material (~80)\*** – a mousse-like viscosity bite material with a high end hardness/rigidity that facilitates trimming; can also be used as a heavy body/tray material.

**Correct VPS No-Flex™†† (~80)\*** – a high viscosity material with a very high durometer specifically made to be used as a tray material in a double arch type impression.

**Correct VPS/Quick Putty (~70)\*** – an easy to mix, non-oily, heavy viscosity putty with a high end hardness to resist distortion and provide maximum hydraulic drive of the light body material.

### Featured Benefits

- **Unsurpassed accuracy** – better fitting restorations
- **Exceptional stability** – molds can be poured multiple times without worrying about distortion
- **Consistent Performance** – various work and set times for your own customized techniques
- **High tear strengths** - easy removal from the mouth
- **No taste or smell** - better patient acceptance
- **Reliability** - proven track record of dependable results for over 10 years

\*Shore A hardness after 1 hour. \*\*Correct VPS Bite Material only. †Tubes available in VPS only, ††Not available in Quick Set.

# Correct VPS® Correct Quick®

Vinyl Polysiloxane  
Impression Material

## Ordering Information

4-Pack cartridge kits include 16 mixing tips; 32-pack bulk cartridge kits include 116 mixing tips.

	Correct VPS®	Correct Quick®
<b>Light Extra Body</b>		
4 – 50ml Cartridge Kit	Q02AC	Q31A
32 – 50ml Cartridge Bulk Pack,	Q09SAC	Q31AB
<i>4-pack also includes 16 Intra-Oral tips; 32-pack also includes 116 Intra-Oral tips</i>		
<b>Light Body</b>		
4 – 50ml Cartridge Kit	Q02A	Q31B
32 – 50ml Cartridge Bulk Pack	Q09SA	Q31BB
1 – 90ml Base/90ml Catalyst Tube Kit	Q01A	N/A
12 – 90ml Base/Catalyst Tube Bulk Pack	Q08A	N/A
<i>4-pack also includes 16 Intra-Oral tips; 32-pack also includes 116 Intra-Oral tips</i>		
<b>Medium Body</b>		
4 – 50ml Cartridge Kit	Q02B	Q31C
32 – 50ml Cartridge Bulk Pack	Q09SB	Q31CB
1 – 90ml Base/Catalyst Tube Kit	Q01B	N/A
12 – 90ml Base/Catalyst Tube Bulk Pack	Q08B	N/A
<b>Heavy Body</b>		
4 – 50ml Cartridge Kit	Q02C	Q31D
32 – 50ml Cartridge Bulk Pack	Q09SC	Q31DB
1 – 90ml Base/Catalyst Tube Kit	Q01C	N/A
12 – 90ml Base/Catalyst Tube Bulk Pack	Q08C	N/A
<b>Correct Bite™ Bite Registration Material</b>		
4 – 50ml Cartridge Kit	Q05	Q31E
32 – 50ml Cartridge Bulk Pack	Q11SA	Q31EB
1 – 90ml Base/Catalyst Tube Kit	Q04	N/A
12 – 90ml Base/Catalyst Tube Bulk Pack	Q10	N/A
<b>No-Flex™ Impression Material</b>		
4 – 50ml Cartridges	Q28AA	N/A
32 – 50ml Cartridge Bulk Pack	Q28SA	N/A
<i>2-pack includes 10 assorted trays and 6 mixing tips</i>		
<b>Putty</b>		
1 Jar 270ml Catalyst & 1 Jar 270ml Base	Q03	Q31G

For more information on our entire line of innovative dental products visit us at [www.pentron.com](http://www.pentron.com) or call us at **1.800.551.0283**

800.551.0283 | 203.265.7397 | [www.pentron.com](http://www.pentron.com)

Correct Plus® | Correct Plus® Fast Set | Correct Quick® | Correct VPS®

Impression

Adhesives

Crown & Bridge

Restoratives

Subject to change or cancellation without notice. Copyright © 2011 Pentron Clinical. All rights reserved.

CQVLT\_Rev01

