



# First-fill®

## Temporary Filling Material

*First-fill* Light Cured Temporary Filling Material is a special polymer based formulation created to address your needs when fabricating small to moderate sized provisional restorations. *First-fill* Light Cured Temporary Filling Material represents an advance over traditional self curing materials due to its semi-flexible nature when set and its ability to be light cured. As a result, you will be able to create provisional restorations that are easy to place, set and remove.



### First-fill® Temporary Filling Material provides:

- > **Convenience** – The single component, syringe dispensed material does not require mixing or the use of matrix bands.
- > **Ease-of-Use** – A putty like consistency makes the product easy to handle, shape and place.
- > **Fast Setting** – First-fill material is a light cured temporary filling material, so it sets upon your command. This eliminates the tedious waiting period required for self curing material.
- > **Bulk Filling** – Because of its translucent color, the product can be light cured in increments up to 5mm.
- > **Easy Removal** – Semi-flexible upon setting enables removal when slight undercuts are present without damaging the margins.

# First-fill®

## Temporary Filling Material

### Indications:

*First-fill* Light Cured Temporary Filling Material is designed for small to moderate sized provisional restorations. It is an ideal temporary filling material for inlays, onlays and endodontic access openings.

### Ordering Information:

#### First Fill Temporary Filling Material Kit

2 - 2ml First Fill Syringes • Instructions and MSDS

REF N19

## Also Try Our Other Great Crown & Bridge Products:

### TempSpan® Provisional System

The *TempSpan* Provisional System is designed to save time by providing a full range of convenient and easy to use products for the fabrication of temporary restorations. *TempSpan* Dual-Cure Temporary Crown & Bridge material is a dual-cure composite available in 7 standard shades that can be light cured intra- and extra-orally. *TempSpan* Clear Matrix material is a transparent vinyl polysiloxane material that enables light activation of the dual cure *TempSpan* C&B material so the temporary is easily fabricated in the matrix. For cementing provisional restorations, *TempSpan* CMT Temporary Cement provides ideal retention, easy clean up and helps reduce sensitivity. Designed to eliminate the time consuming step of polishing provisional restorations, *TempSpan* Glaze creates a lustrous gloss on provisional restorations and is compatible with both halogen and LED curing lights.



For more information on our entire line of innovative dental products visit us at [www.pentron.com](http://www.pentron.com) or call us at **1.800.551.0283**

800.551.0283 | [www.pentron.com](http://www.pentron.com) | 203.265.7397

TempSpan® | Build-It® FR™ | FibreKleer® Post | First-Fill® Crown & Bridge Impressions Adhesives Restoratives

Subject to change or cancellation without notice. Copyright © 2009 Pentron Clinical Technologies, LLC. All rights reserved.

FLT\_Rev01

