



Bond-1®/Bond-1® C&B™ Primer/Adhesive

Bond-1 Primer/Adhesive, Pentron's single-bottle moist dentin bonding system, is easy to use, provides high bond strengths and helps eliminate sensitivity. Proven to create an ideal sealed surface and interpenetrated network of the polymerized primers and unfilled resin, *Bond-1* Primer/Adhesive yields high bond strengths with film thickness of only 8µm. Its ease of application, combined with its low film thickness, assures the ideal surface to bond to all types of restorative materials. *Bond-1* Primer/Adhesive is available in both an acetone and alcohol-based formula to accommodate all of your bonding needs.

Features & Benefits

- 5th generation, moist dentin bonding system
- High bond strengths
- Bonds to virtually all direct & indirect restorations
- Dual-cure activator allows for compatibility with most self- and dual-cure resin materials
- Seals dentin completely by forming a hybrid zone for optimal adhesion and bonding
- Helps eliminate post-operative sensitivity
- Film thickness of only **8 microns**
- Choice of solvent base:
 - Alcohol provides longer working time
 - Acetone dries quicker, saving time
- Save up to **50%** over the leading 5th generation bonding agent brands

Bond-1® [acetone-based formula]



Bond-1® C&B™ [alcohol-based formula]

Bond-1[®] Bond-1[®] C&B[™]

Primer/Adhesive

Ordering Information:

Bond-1[®] Kit REF N011

- 1 – 4ml bottle *Bond-1* Primer/Adhesive
- 1 – 3ml bottle *Bond-1* Dual Cure Activator
- 2 – 5ml syringe 37% Etch gel
- 1 – Mixing Well
- 20 – 22 gauge needle tips
- 50 – Disposable Applicators
- Instructions

Bond-1 Refills

- | | |
|--|------------|
| 1 – 4ml bottle <i>Bond-1</i> Primer/Adhesive | REF N011AA |
| 1 – 6ml bottle <i>Bond-1</i> Primer/Adhesive | REF N011AB |
| 1 – 3ml bottle <i>Bond-1</i> Dual Cure Activator | REF N011AC |

Bond-1 C&B Refills

- | | |
|--|----------|
| 1 – 6ml bottle <i>Bond-1</i> C&B Primer/Adhesive | REF N31B |
|--|----------|

Additional Kits and Refills:

- | | |
|--|-----------|
| 2 – 5ml syringe 37% Etch gel w/20-22 g. tips | REF N011B |
|--|-----------|

For more information on our entire line of innovative dental products visit us at www.pentron.com or call us at **1.800.551.0283**

800.551.0283 | 203.265.7397 | www.pentron.com