

Correct Plus®

Hydrophilic
Impression Materials

Ordering Information

Correct Plus Impression Material Cartridges and Putty

- 4 - 50ml Cartridges, Instructions and Accessories
- 32 - 50ml Cartridges, Instructions and Accessories
- 2 - 380ml Mighty Mix Cartridges, 20 Mixing Tips, 2 Bayonet Rings and Instructions
- 1 - 230ml Jar each - Putty Base & Catalyst, Instructions and Accessories

Available in the following viscosities:

Unflavored	4 Cartridges	32 Cartridges	Mighty Mix
Light Body	REF Q34BA	REF Q34BC	—
Universal™ Body	REF Q34DA	REF Q34DC	REF Q34MA
Thick'n Thin® Light Body	REF Q34GA	REF Q34GC	—
Auto-Mix Putty	REF Q34HA	REF Q34HC	REF Q34MB
Bite Superfast Material	REF Q34IA	REF Q34IC	—
Tray Material	REF Q34JA	REF Q34JC	REF Q34MC
Putty — Base (230ml) & Catalyst (230ml)		REF Q34H	

Flavored	4 Cartridges	32 Cartridges	
Berry Tray Material	REF Q34JAA	REF Q34JAC	
Berry Auto-Mix Putty	REF Q34HAA	REF Q34HAC	
Berry Putty - Base (230ml) & Catalyst (230ml)		REF Q34HCA	

Correct Plus Fast Set Impression Material Cartridges

- 4 - 50ml Cartridges, Instructions and Accessories
- 32 - 50ml Cartridges, Instructions and Accessories
- 2 - 380ml *Mighty Mix* Cartridges, 20 Mixing Tips, 2 Bayonet Rings and Instructions

Available in the following viscosities:

Unflavored	4 Cartridges	32 Cartridges	Mighty Mix
Thick n' Thin® Light Body	REF Q36C	REF Q36CC	
Universal™ Body	REF Q36B	REF Q36BB	REF Q36MA
Auto-Mix Putty	REF Q36A	REF Q36AA	REF Q36MB

Impression Mixing Tips & Accessories

100 - Refill Tips (Green)	REF Q19BB	50 - Dynamic Mixing Tips	REF Q34MM
100 - Refill Tips (Purple)	REF Q19CB	2 - Bayonet Rings	REF Q34MR
100 - Refill Tips (Pink)	REF Q19AB		
100 - Intra-oral Tips	REF Q16		
50ml - Dispenser Gun	REF Q19		

For more information on our entire line of innovative dental products visit us at www.pentron.com or call us at 1.800.551.0283

Correct Plus®

Hydrophilic
Impression Materials



Sybron Dental Specialties

53 North Plains Industrial Road
Wallingford, CT 06492

t | 800.551.0283 | 203.265.7397
w | www.pentron.com



Correct Plus® Hydrophilic Impression Materials

The road to a perfect impression starts here - you choose how long it takes to get there.

This extensive line of hydrophilic impression material is available in a variety of viscosities with 2 set times: Standard and Fast Set. *Correct Plus* Impression Materials are designed to consistently provide highly accurate and detailed impressions every time you use them. Its versatile setting times, multiple viscosities, and choice of 50ml auto-mix cartridges or 380ml Mighty Mix™ cartridges designed for automated dispensing machines make *Correct Plus* Impression Material the perfect choice for fabricating crowns, bridges, inlays, onlays, veneers, and more.



- > Low contact angle of $\approx 30^\circ$ for accurate, reliable impressions in a wet or dry environment
- > Choice of work/set times for single unit to full arch impressions
- > Outstanding tear strength ensures reliable impressions
- > Minimized outgassing time for immediate pour of models
- > Unique viscosities including Thick n' Thin® Light Body and Auto-Mix Putty
- > Choice of 50ml auto-mix cartridges or 380ml *Mighty-Mix* Cartridges
- > Reduced hand fatigue and time to fill the tray



Thick n' Thin Light Body initially stacks, then melts into the sulcus.

Correct Plus® Hydrophilic Impression Material

Correct Plus Hydrophilic Impression Material is a versatile impression material designed to accommodate all impression techniques. A low contact angle of $\approx 30^\circ$ guarantees outstanding detail in the presence of fluids. Exceptional tear strength ensures your detailed impression will remain intact upon removal. With a working time of 1 minute and 30 seconds and an oral setting time of 3 minutes and 30 seconds, *Correct Plus* is the ideal material for full arch impressions or any situation where a longer set time is desired. For optimal patient acceptance, *Correct Plus* Impression Material is available in unflavored and Berry.



Correct Plus Standard
1:30 Working Time
3:30 Intraoral Set Time

Correct Plus® Hydrophilic Impression Material FastSet

Correct Plus Fast Set
1:00 Working Time
2:00 Intraoral Set Time

Formulated to set in the mouth in 2 minutes, *Correct Plus* Fast Set Impression Material is the perfect choice for exceptionally fast and accurate impressions of 1 – 3 units. Formulated with a unique nano-filler, *Correct Plus* Fast Set Impression Material yields a tear strength that is 11% higher than standard *Correct Plus* Impression Material. This new formulation improves upon the already impressive physical properties of standard *Correct Plus* Impression Material further minimizing the chance of tearing and ensuring that impressions are dimensionally stable pour after pour. A low contact angle of $\approx 30^\circ$ minimizes the occurrence of voids and ensures accuracy, especially in a wet environment. *Correct Plus* Fast Set Impression Material is available in our three most popular viscosities, *Thick n' Thin* Light Body, *Universal*™ Body, and *Auto-Mix Putty*, to accurately support three-in-one tray, stock tray, and custom tray impression techniques.



Correct Plus Mighty Mix Cartridge System for use in automated mixing machines like Pentamix®†

"Fast and easy dispensing with better economics than when using cartridges to fill the tray"

† Pentamix is a registered trademark of 3M ESPE.

- **Gregori M. Kurtzman, DDS, MAGD, FACD, DICOI**

Correct Plus Hydrophilic Impression Material Selection Chart

Viscosity	Description	Indication*	Set Times	Color	Delivery Options
<i>Thick n' Thin</i> Light Body	A unique, thixotropic light viscosity material that initially stacks on the tooth when syringed then flows into the sulcus to capture critical subgingival detail.	W	Standard & Fast Set	Green	1:1
Light Body	A traditional light viscosity wash material with exceptional flowability.	W	Standard	Purple	1:1
<i>Universal</i> Body	Truly an all-purpose impression material. Use it as a tray material and/or a light body. A perfect substitute for alginate for night guards, bleach trays, or any other impression procedure.	W & C	Standard & Fast Set	Blue	1:1, 5:1
Auto-Mix Putty	A heavy viscosity material that provides the hydraulic drive of a hand-mix putty in a convenient auto-mix cartridge.	Q, FA, TIO	Standard & Fast Set	Cyan	1:1, 5:1
Tray Material	A heavy viscosity material with a mousse-like consistency and high end hardness.	Q, FA, TIO	Standard	Dark Purple	1:1, 5:1
Putty	An easy to mix, heavy viscosity putty with a high end hardness to resist distortion and provide maximum hydraulic drive of the light body material.	Q & FA	Standard	Light Blue	Hand Mix
Bite Superfast Material	A superfast setting bite registration material with a mousse-like viscosity. When set, its rigidity facilitates trimming.	Bite Registration	Fast Set	Dark Blue	1:1

*Indication Abbreviations: W=Wash, C=Custom Tray, Q=Quadrant Tray, FA=Full Arch Tray, TIO=Three-In-One Impression Tray
1:1=50ml 1:1 Auto-Mix Cartridges / 5:1=380ml 5:1 *Mighty Mix* Cartridges