



PrepEze®

Desensitizer

PrepEze Desensitizer is designed to eliminate post-operative sensitivity! This colorless solution can be applied in just 30 seconds, does not require separate curing, and is non-irritating to soft tissue. *PrepEze* Desensitizer is compatible with all conventional cements, amalgam and bonding agents.

Effectively eliminate sensitivity in all areas of dentistry!

PrepEze Desensitizer is indicated for the treatment of exposed cervical dentin, post-operative desensitization, placement on crown preps prior to temporization, moistening the dentin surface after etching and use with scaling and root planing.

PrepEze Desensitizer eliminates sensitivity via polymerization in combination with the bonding agent in the dentin tubules, creating a physiological seal. *PrepEze* Desensitizer contains benzalkonium chloride, a fast-acting anti-microbial agent to rid the dentinal surface of residual bacteria (a major cause of sensitivity), HEMA, a well-proven wetting agent, and sodium fluoride, an established desensitizer.

Featured Benefits

- Easy to apply
- Non irritating to soft tissues
- Compatible with all conventional cements, amalgam and bonding agents
- Eliminates sensitivity by creating a physiological seal in dentin tubules
- Contains benzalkonium chloride, an anti-microbial agent to eliminate residual bacteria
- Contains fluoride
- Ideal for the treatment of exposed cervical dentin
- Indicated for routine post-operative use
- Ideal for treating crown preps prior to temporization
- Moistening the dentin surface after etching
- Eliminates sensitivity after scaling and root planing

PrepEze®

Desensitizer

Ordering Information:

1 - 5ml bottle PrepEze® Desensitizer

REF N33B

For sensitivity-free cementation, try with:

Mojo™ Veneer Cement Kit

REF N06K

8 - 1.6gm syringes - Cement and Try-in,
1 each shade: Clear, Dark, Light and Movie Star White
1 - 4ml bottle *Bond-1* Primer/ Adhesive
1 - 1ml Syringe 37% Phosphoric Acid Etching Gel
1 - 3ml bottle Silane
40 Plastic tips, Instructions

Lute-It® with Bond-1® Kit

REF N04

8 – 3.3gm syringes Luting Cement - Base:
Light, Dark, Clear, Universal, Grey, Opac. White, Opac. Universal, Dark Yellow
8 – 3ml syringes Try-in Gel:
Light, Dark, Clear, Universal, Grey, Opac. White, Opac. Universal, Dark Yellow
2 – 3.3gm syringes Catalyst
1 – 1ml syringe 37% Etch gel
1 – 3ml bottle Silane
1 – 6ml bottle *Bond-1* Primer/Adhesive
Instructions and Accessories

Lute-It® with Bond-1® Intro Kit

REF N04B

3 – 3.3gm syringes Luting Cement - Base:
Light Yellow, Clear, Universal
1 – 1ml syringe Stain:
Opaque Universal
1 – 4ml bottle *Bond-1* Primer/Adhesive
1 – 3ml bottle Silane
1 – 3.3gm syringe Catalyst
1 – 1ml syringe 37% Etch gel
Instructions and Accessories

Cement-It® Universal C&B™ Resin Cementation System

REF N33

1 - 4ml Auto-mix Syringe A2 catalyst & base
1 - 4ml Auto-mix Syringe Opaque White catalyst & base
1 - 6ml bottle *Bond-1* C&B Primer/Adhesive
1 - 5ml bottle *PrepEze* Desensitizer
1 - 3ml syringe of 37% Etching gel
1 - 3ml bottle of Silane
Accessories and Instructions

For more information on our entire line of innovative dental products and refills, visit us at www.pentron.com or call us at **1.800.551.0283**

800.551.0283 | 203.265.7397 | www.pentron.com

Artiste® | Simile® | Flow-It® ALC™ | Alert® | Fini™ | Fusio™ | PrepEze®

Restoratives

Impression

Crown & Bridge

Adhesives

* Subject to change or cancellation without notice. Copyright © 2010 Pentron Clinical. All rights reserved.

SLT_Rev02

 **Pentron
Clinical**
Sybron Dental Specialties