



Protect-It[®]

Composite Surface Sealant



Protect-It Composite Surface Sealant is a protective resin used to seal microcracks and surface irregularities present in composite restorations after finishing and polishing. The protective resin surface helps decrease composite wear so your restorations last longer and seals margins to help resist microleakage.

Protect-It[™] Featured benefits:

- Protects against microleakage
- Helps decrease the wear of posterior composites
- Penetrates into polished composite
- Easy maintenance – can be reapplied on follow-up visits
- Simple application that requires no refinishing or repolishing



Fini[™]

Polishing Disk System



Fini Polishing Disks provide an incredibly smooth and brilliant polish to all composites. The innovative centerless disk provides greater access and increases the usable surface area of each disk, making finishing and polishing procedures quick and easy. Selected disks are designed with grit on both sides requiring no disk change for lingual access.

Fini[™] Polishing Disk Featured benefits:

- Centerless design eliminates gouging and discoloration common in other disk systems
- Disks are color coded by grit range for easy identification
- Disks are thin and flex easy for very good interproximal and marginal access

Protect-It[®]

Composite Surface Sealant

Ordering Information

Protect-It Composite Surface Sealant Kit **N10**

Contents:

- 1 - 7ml bottle *Protect-It* Composite Surface Sealant
- 2 - 5ml syringes 37% Etch Gel
- Instructions and Accessories

Refills:

- 1 - 7ml bottle *Protect-It* Composite Surface Sealant N10A
- 2 - 5ml syringes 37% Etch Gel N01IB
- 100 - Disposable Applicators N10B
- 50 - 22 gauge bent needle tips T26
- 100- 22 gauge bent needle tips T26B

Ordering Information

Fini Polishing Disks Kit **N34**

Kit Contains:

- 180 disks (20 ea. of nine sizes)
- Instructions and Accessories

Refills:

- 50 - Black Course - Large Outer Grit N34AA
- 50 - Black Course - Small Outer Grit N34AB
- 50 - Blue Medium - Large Outer Grit N34AC
- 50 - Blue Medium - Large Inner Grit N34AD
- 50 - Blue Medium - Small Outer Grit N34AE
- 50 - Red Fine - Large Grit both Sides N34AF
- 50 - Red Fine - Small Grit both Sides N34AG
- 50 - White Extra Fine - Large Grit both Sides N34AH
- 50 - White Extra Fine - Small Grit both Sides N34AI
- 2 - Contra Angle Mandrels N34AJ

For more information on our entire line of innovative dental products visit us at www.pentron.com or call us at **1.800.551.0283**

Fini[™]

Polishing Disks



800.551.0283 | 203.265.7397 | www.pentron.com

Artiste[®] | Simile[®] | Flow-It[®]ALC[™] | Alert[®] | Fini[™] | Protect- It[®]

Restoratives

Impression

Crown & Bridge

Adhesives



Sybron Dental Specialties

* Subject to change or cancellation without notice. Copyright © 2010 Pentron Clinical. All rights reserved.